



LIF Antibody

Product Code	CSB-PA012928ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15018
Immunogen	Recombinant Human Leukemia inhibitory factor protein (23-202AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukemia inhibitory factor (LIF) (Differentiation-stimulating factor) (D factor) (Melanoma-derived LPL inhibitor) (MLPLI) (Emfilermin), LIF, HILDA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LIF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.