




ULBP2 Antibody

Product Code	CSB-PA880151ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZM5
Immunogen	Recombinant Human NKG2D ligand 2 protein (26-216AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	UL16-binding protein 2 (ALCAN-alpha) (NKG2D ligand 2) (N2DL-2) (NKG2DL2) (Retinoic acid early transcript 1H), ULBP2, N2DL2 RAET1H
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ULBP2

Image



Immunohistochemistry of paraffin-embedded human duodenum tissue using CSB-PA880151ESR1HU at dilution of 1:100