




TACSTD2 Antibody

Product Code	CSB-PA023072ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09758
Immunogen	Recombinant Human Tumor-associated calcium signal transducer 2 protein (27-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor-associated calcium signal transducer 2 (Cell surface glycoprotein Trop-2) (Membrane component chromosome 1 surface marker 1) (Pancreatic carcinoma marker protein GA733-1), TACSTD2, GA733-1 M1S1 TROP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	TACSTD2

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA023072ESR1HU at dilution of 1:100