




PKD2L1 Antibody

Product Code	CSB-PA865199DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0L9
Immunogen	Recombinant Human Polycystic kidney disease 2-like 1 protein (561-805AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycystic kidney disease 2-like 1 protein (Polycystin-2 homolog) (Polycystin-2L1) (Polycystin-L) (Polycystin-L1), PKD2L1, PKD2L PKDL TRPP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PKD2L1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA865199DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.