




PRKAB2 Antibody

Product Code	CSB-PA527326ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43741
Immunogen	Recombinant Human 5'-AMP-activated protein kinase subunit beta-2 protein (1-272AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	5'-AMP-activated protein kinase subunit beta-2 (AMPK subunit beta-2), PRKAB2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PRKAB2

Image



Western blot

All lanes: PRKAB2 antibody at 5µg/ml

Lane 1: Rat brain tissue

Lane 2: Mouse heart tissue

Lane 3: Mouse gonad tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 31, 22 kDa

Observed band size: 31 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.