




# NOV Antibody

<b>Product Code</b>	CSB-PA015956ESR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P48745
<b>Immunogen</b>	Recombinant Human Protein NOV homolog protein (32-357AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein NOV homolog (NovH) (CCN family member 3) (Insulin-like growth factor-binding protein 9) (IBP-9) (IGF-binding protein 9) (IGFBP-9) (Nephroblastoma-overexpressed gene protein homolog), NOV, CCN3 IGFBP9 NOVH
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	NOV

## Image



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA015956ESR2HU at dilution of 1:100