




NOV Antibody

Product Code	CSB-PA015956ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48745
Immunogen	Recombinant Human Protein NOV homolog protein (32-357AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NOV homolog (NovH) (CCN family member 3) (Insulin-like growth factor-binding protein 9) (IBP-9) (IGF-binding protein 9) (IGFBP-9) (Nephroblastoma-overexpressed gene protein homolog), NOV, CCN3 IGFBP9 NOVH
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NOV

Image



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA015956ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.