




PRKAG3 Antibody

Product Code	CSB-PA883379ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGI9
Immunogen	Recombinant Human 5'-AMP-activated protein kinase subunit gamma-3 protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	5'-AMP-activated protein kinase subunit gamma-3 (AMPK gamma3) (AMPK subunit gamma-3), PRKAG3, AMPKG3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKAG3

Image



Western blot

All lanes: Carbonic anhydrase 1 antibody at 6µg/ml + Rat brain tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 55, 52 kDa

Observed band size: 55 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.