




# ADAP1 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA001322DSR1HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | O75689   |
| <b>Immunogen</b>           | Recombinant Human Arf-GAP with dual PH domain-containing protein 1 protein (1-374AA)   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA, IHC; Recommended dilution: IHC:1:20-1:200   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.  |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Arf-GAP with dual PH domain-containing protein 1 (Centaurin-alpha-1) (Cnt-a1) (Putative MAPK-activating protein PM25), ADAP1, CENTA1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Neuroscience   |
| <b>Target Names</b>        | ADAP1  |

## Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA001322DSR1HU at dilution of 1:100