





CDK10 Antibody

Product Code	CSB-PA005060ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15131
Immunogen	Recombinant Human Cyclin-dependent kinase 10 protein (1-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 10 (EC 2.7.11.22) (Cell division protein kinase 10) (Serine/threonine-protein kinase PISLRE), CDK10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDK10

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA005060ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA005060ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.