



NDUFS7 Antibody

Product Code	CSB-PA015666ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75251
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial protein (39-213AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial (EC 1.6.5.3) (EC 1.6.99.3) (Complex I-20kD) (CI-20kD) (NADH-ubiquinone oxidoreductase 20 kDa subunit) (PSST subunit), NDUFS7
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NDUFS7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.