




AHSP Antibody

Product Code	CSB-PA878935ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZD4
Immunogen	Recombinant Human Alpha-hemoglobin-stabilizing protein (1-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-hemoglobin-stabilizing protein (Erythroid differentiation-related factor) (Erythroid-associated factor), AHSP, EDRF ERAF
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	AHSP

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA878935ESR2HU at dilution of 1:100