





ORC4 Antibody

Product Code	CSB-PA017234ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43929
Immunogen	Recombinant Human Origin recognition complex subunit 4 protein (187-436AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Origin recognition complex subunit 4, ORC4, ORC4L
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ORC4

Image



Western blot
 All lanes: ORC4 antibody at 4µg/ml
 Lane 1: HeLa whole cell lysate
 Lane 2: Jurkat whole cell lysate
 Lane 3: A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 51, 41, 43 kDa
 Observed band size: 51 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA017234ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.