





DEFB4A Antibody

Product Code	CSB-PA006705ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15263
Immunogen	Recombinant Human Beta-defensin 4A protein (1-64AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-defensin 4A (Beta-defensin 2) (BD-2) (hBD-2) (Defensin, beta 2) (Skin-antimicrobial peptide 1) (SAP1), DEFB4A; DEFB4B, DEFB102 DEFB2 DEFB4;
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	DEFB4A

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA006705ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA006705ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.