





DRD3 Antibody

Product Code	CSB-PA007180ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35462
Immunogen	Recombinant Human D(3) dopamine receptor protein (213-329AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	D(3) dopamine receptor (Dopamine D3 receptor), DRD3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DRD3

Image



Western blot
 All lanes: DRD3 antibody at 6µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 45, 41 kDa
 Observed band size: 45 kDa



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA007180ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.