




HYLS1 Antibody

Product Code	CSB-PA846663DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96M11
Immunogen	Recombinant Human Hydrolethalus syndrome protein 1 protein (70-299AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydrolethalus syndrome protein 1, HYLS1, HLS
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HYLS1

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA846663DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.