





COIL Antibody

Product Code	CSB-PA005713DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38432
Immunogen	Recombinant Human Coilin protein (297-576AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coilin (p80-coilin), COIL, CLN80
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	COIL

Image



Western blot
 All lanes: COIL antibody at 3µg/ml
 Lane 1: MCF-7 whole cell lysate
 Lane 2: HeLa whole cell lysate
 Lane 3: HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 63 kDa
 Observed band size: 63 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA005713DSR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.