



groL Antibody, HRP conjugated

Product Code	CSB-PA394210HB01YBL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4TRR0
Immunogen	Recombinant Yersinia pestis 60 kDa chaperonin protein (1-548AA)
Raised In	Rabbit
Species Reactivity	Yersinia pestis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60 kDa chaperonin (GroEL protein) (Protein Cpn60), groL, groEL
Immunogen Species	Yersinia pestis
Research Area	Others
Target Names	groL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.