




STX11 Antibody

Product Code	CSB-PA022889ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75558
Immunogen	Recombinant Human Syntaxin-11 protein (1-287AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-11, STX11
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	STX11

Image



Western blot


All lanes: STX11 antibody at 6µg/ml + Mouse skeletal muscle tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 33 kDa

Observed band size: 33 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA022889ESR2HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.