




CAPZB Antibody

Product Code	CSB-PA004514ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47756
Immunogen	Recombinant Human F-actin-capping protein subunit beta protein (1-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-actin-capping protein subunit beta (CapZ beta), CAPZB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAPZB

Image



Western blot

All lanes: CAPZB antibody at 6µg/ml

Lane 1: 293T whole cell lysate

Lane 2: A431 whole cell lysate


Lane 3: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 32, 31 kDa

Observed band size: 32, 31 kDa



Immunohistochemistry of paraffin-embedded

human testis tissue using CSB-

PA004514ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.