




SHC4 Antibody

Product Code	CSB-PA747568ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6S5L8
Immunogen	Recombinant Human SHC-transforming protein 4 protein (1-185AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	SHC-transforming protein 4 (Rai-like protein) (RaLP) (SHC-transforming protein D) (hShcD) (Src homology 2 domain-containing-transforming protein C4) (SH2 domain protein C4), SHC4, SHCD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SHC4

Image



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA747568ESR2HU at dilution of 1:100