




STXBP3 Antibody

Product Code	CSB-PA022907ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00186
Immunogen	Recombinant Human Syntaxin-binding protein 3 protein (343-592AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-binding protein 3 (Platelet Sec1 protein) (PSP) (Protein unc-18 homolog 3) (Unc18-3) (Protein unc-18 homolog C) (Unc-18C), STXBP3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	STXBP3

Image



Western blot

All lanes: STXBP3 antibody at 10µg/ml

Lane 1: HepG2 whole cell lysate


Lane 2: K562 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 68 kDa

Observed band size: 68 kDa



Immunohistochemistry of paraffin-embedded

human pancreatic tissue using CSB-

PA022907ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.