




AP1M2 Antibody

Product Code	CSB-PA895771ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6Q5
Immunogen	Recombinant Human AP-1 complex subunit mu-2 protein (164-423AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	AP-1 complex subunit mu-2 (AP-mu chain family member mu1B) (Adaptor protein complex AP-1 subunit mu-2) (Adaptor-related protein complex 1 subunit mu-2) (Clathrin assembly protein complex 1 mu-2 medium chain 2) (Golgi adaptor HA1/AP1 adaptin mu-2 subunit) (Mu-adaptin 2) (Mu1B-adaptin), AP1M2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AP1M2

Image



Western blot

All lanes: AP1M2 antibody at 2μg/ml


Lane 1: Mouse small intestine tissue

Lane 2: Mouse kidney tissue


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48 kDa

Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA895771ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA895771ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.