




MAPK4 Antibody

Product Code	CSB-PA013457ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31152
Immunogen	Recombinant Human Mitogen-activated protein kinase 4 protein (348-587AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 4 (MAP kinase 4) (MAPK 4) (EC 2.7.11.24) (Extracellular signal-regulated kinase 4) (ERK-4) (MAP kinase isoform p63) (p63-MAPK), MAPK4, ERK4 PRKM4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAPK4

Image



Western blot

All lanes: MAPK4 ntibody at 6µg/ml

Lane 1: Mouse brain tissue

Lane 2: 293T whole cell lysate


Lane 3: Jurkat whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 66 kDa

Observed band size: 66 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA013457ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA013457ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.