




# RPH3AL Antibody

<b>Product Code</b>	CSB-PA891975DSR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UNE2
<b>Immunogen</b>	Recombinant Human Rab effector Noc2 protein (161-315AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Rab effector Noc2 (No C2 domains protein) (Rabphilin-3A-like protein), RPH3AL, NOC2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	RPH3AL

## Image



### Western blot

All lanes: RPH3AL antibody at 12µg/ml

Lane 1: LO2 whole cell lysate


Lane 2: K562 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35, 32 kDa

Observed band size: 35 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA891975DSR1HU at dilution of 1:100



---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.