




RPH3AL Antibody

Product Code	CSB-PA891975DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNE2
Immunogen	Recombinant Human Rab effector Noc2 protein (161-315AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab effector Noc2 (No C2 domains protein) (Rabphilin-3A-like protein), RPH3AL, NOC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RPH3AL

Image



Western blot

All lanes: RPH3AL antibody at 12μg/ml

Lane 1: LO2 whole cell lysate


Lane 2: K562 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35, 32 kDa

Observed band size: 35 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA891975DSR1HU at dilution of 1:100

