




IL27 Antibody

Product Code	CSB-PA847702ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEV9
Immunogen	Recombinant Human Interleukin-27 subunit alpha protein (29-243AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-27 subunit alpha (IL-27 subunit alpha) (IL-27-A) (IL27-A) (Interleukin-30) (p28), IL27, IL27A IL30
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL27

Image



Western blot
 All lanes: IL27 antibody at 2µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 27 kDa
 Observed band size: 27 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.