



RGS3 Antibody

Product Code	CSB-PA019655ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49796
Immunogen	Recombinant Human Regulator of G-protein signaling 3 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of G-protein signaling 3 (RGP3) (RGS3), RGS3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.