



SNAP91 Antibody

Product Code	CSB-PA021876DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60641
Immunogen	Recombinant Human Clathrin coat assembly protein AP180 protein (628-907AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Clathrin coat assembly protein AP180 (91 kDa synaptosomal-associated protein) (Clathrin coat-associated protein AP180) (Phosphoprotein F1-20), SNAP91, KIAA0656
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNAP91
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.