




# FABP5 Antibody

<b>Product Code</b>	CSB-PA007946ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q01469
<b>Immunogen</b>	Recombinant Human Fatty acid-binding protein, epidermal protein (1-135AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Fatty acid-binding protein, epidermal (Epidermal-type fatty acid-binding protein) (E-FABP) (Fatty acid-binding protein 5) (Psoriasis-associated fatty acid-binding protein homolog) (PA-FABP), FABP5
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	FABP5

## Image



Western blot

All lanes: FABP5 antibody at 5µg/ml + HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 15 kDa

Observed band size: 15 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.