





CCL11 Antibody

Product Code	CSB-PA004775ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51671
Immunogen	Recombinant Human Eotaxin protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eotaxin (C-C motif chemokine 11) (Eosinophil chemotactic protein) (Small-inducible cytokine A11), CCL11, SCYA11
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCL11

Image



Western blot
 All lanes: CCL11 antibody at 5µg/ml + Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 11 kDa
 Observed band size: 11 kDa



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA004775ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.