



CNTNAP2 Antibody

Product Code	CSB-PA887030ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHC6
Immunogen	Recombinant Human Contactin-associated protein-like 2 protein (180-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Contactin-associated protein-like 2 (Cell recognition molecule Caspr2), CNTNAP2, CASPR2 KIAA0868
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CNTNAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.