





SLC4A2 Antibody

Product Code	CSB-PA021667DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04920
Immunogen	Recombinant Human Anion exchange protein 2 protein (95-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anion exchange protein 2 (AE 2) (Anion exchanger 2) (Non-erythroid band 3-like protein) (BND3L) (Solute carrier family 4 member 2), SLC4A2, AE2 EPB3L1 HKB3 MPB3L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC4A2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA021667DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA021667DSR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.