




# TERF1 Antibody

|                            |                                                                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA023386ESR2HU                                                                                                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                             |
| <b>Uniprot No.</b>         | P54274                                                                                                                                                    |
| <b>Immunogen</b>           | Recombinant Human Telomeric repeat-binding factor 1 protein (260-439AA)                                                                                   |
| <b>Raised In</b>           | Rabbit                                                                                                                                                    |
| <b>Species Reactivity</b>  | Human                                                                                                                                                     |
| <b>Tested Applications</b> | ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200                                                                                     |
| <b>Form</b>                | Liquid                                                                                                                                                    |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                                            |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                                                                                         |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                                                                                 |
| <b>Isotype</b>             | IgG                                                                                                                                                       |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                |
| <b>Alias</b>               | Telomeric repeat-binding factor 1 (NIMA-interacting protein 2) (TTAGGG repeat-binding factor 1) (Telomeric protein Pin2/TRF1), TERF1, PIN2 TRBF1 TRF TRF1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                      |
| <b>Research Area</b>       | Epigenetics and Nuclear Signaling                                                                                                                         |
| <b>Target Names</b>        | TERF1                                                                                                                                                     |

## Image



### Western blot

All lanes: TERF1 antibody at 2μg/ml

Lane 1: k562 whole cell lysate

Lane 2: 293T whole cell lysate


Lane 3: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 51, 49 kDa

Observed band size: 51 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA023386ESR2HU at dilution of 1:100

---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.