



XCL1 Antibody

Product Code	CSB-PA026186ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47992
Immunogen	Recombinant Human Lymphotactin protein (22-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphotactin (ATAC) (C motif chemokine 1) (Cytokine SCM-1) (Lymphotaxin) (SCM-1-alpha) (Small-inducible cytokine C1) (XC chemokine ligand 1), XCL1, LTN SCYC1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	XCL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.