




STXBP4 Antibody

Product Code	CSB-PA757971ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZWJ1
Immunogen	Recombinant Human Syntaxin-binding protein 4 protein (344-553AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-binding protein 4 (Syntaxin 4-interacting protein) (STX4-interacting protein) (Synip), STXBP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STXBP4

Image



Western blot

All lanes: STXBP4 antibody at 10µg/ml + PC-3 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 62, 27 kDa
 Observed band size: 62 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.