





ACRV1 Antibody

Product Code	CSB-PA001186DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26436
Immunogen	Recombinant Human Acrosomal protein SP-10 protein (22-265AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acrosomal protein SP-10 (Acrosomal vesicle protein 1), ACRV1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACRV1

Image



Western blot
 All lanes: ACRV1 antibody at 4μg/ml
 Lane 1: Mouse liver tissue
 Lane 2: Mouse gonad tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 29, 27, 23, 21, 20, 19, 18, 17, 14, 9 kDa
 Observed band size: 29 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA001186DSR1HU at dilution of 1:100