



MAVS Antibody

Product Code	CSB-PA772018ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z434
Immunogen	Recombinant Human Mitochondrial antiviral-signaling protein (1-65AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial antiviral-signaling protein (MAVS) (CARD adapter inducing interferon beta) (Cardif) (Interferon beta promoter stimulator protein 1) (IPS-1) (Putative NF-kappa-B-activating protein 031N) (Virus-induced-signaling adapter) (VISA), MAVS, IPS1 KIAA1271 VISA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAVS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.