





APOBEC3C Antibody

Product Code	CSB-PA889097ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRW3
Immunogen	Recombinant Human DNA dC->dU-editing enzyme APOBEC-3C protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA dC->dU-editing enzyme APOBEC-3C (A3C) (EC 3.5.4.-) (APOBEC1-like) (Phorbolin I), APOBEC3C, APOBEC1L PBI
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	APOBEC3C

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA889097ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA889097ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.