




ATRIP Antibody

Product Code	CSB-PA848837ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXE1
Immunogen	Recombinant Human ATR-interacting protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATR-interacting protein (ATM and Rad3-related-interacting protein), ATRIP, AGS1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ATRIP

Image



Western blot
 All lanes: ATRIP antibody at 2µg/ml + PC-3 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 86, 84, 77, 73 kDa
 Observed band size: 86 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.