




CA8 Antibody

Product Code	CSB-PA004379ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35219
Immunogen	Recombinant Human Carbonic anhydrase-related protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbonic anhydrase-related protein (CARP) (Carbonic anhydrase VIII) (CA-VIII), CA8, CALS
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	CA8

Image



Western blot

All lanes: CA8 antibody at 4µg/ml

Lane 1: Mouse skeletal muscle tissue


Lane 2: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 33 kDa

Observed band size 33 kDa



Immunohistochemistry of paraffin-embedded

human breast cancer using CSB-

PA004379ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.