




NCAPH Antibody

Product Code	CSB-PA621857DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15003
Immunogen	Recombinant Human Condensin complex subunit 2 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Condensin complex subunit 2 (Barren homolog protein 1) (Chromosome-associated protein H) (hCAP-H) (Non-SMC condensin I complex subunit H) (XCAP-H homolog), NCAPH, BRRN BRRN1 CAPH KIAA0074
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NCAPH

Image


Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA621857DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.