




SATB1 Antibody

Product Code	CSB-PA020719DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01826
Immunogen	Recombinant Human DNA-binding protein SATB1 protein (484-763AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein SATB1 (Special AT-rich sequence-binding protein 1), SATB1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SATB1

Image



Western blot

All lanes: SATB1 antibody at 2µg/ml

Lane 1: Mouse brain tissue

Lane 2: 293T whole cell lysate

Lane 3: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 86, 90 kDa

Observed band size: 86 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.