




NTF3 Antibody

Product Code	CSB-PA016120ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20783
Immunogen	Recombinant Human Neurotrophin-3 protein (19-257AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurotrophin-3 (NT-3) (HDNF) (Nerve growth factor 2) (NGF-2) (Neurotrophic factor), NTF3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NTF3

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA016120ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.