



DAPK1 Antibody

Product Code	CSB-PA006499ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53355
Immunogen	Recombinant Human Death-associated protein kinase 1 protein (1141-1430AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Death-associated protein kinase 1 (DAP kinase 1) (EC 2.7.11.1), DAPK1, DAPK
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DAPK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.