




TRPM2 Antibody

Product Code	CSB-PA025084ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94759
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily M member 2 protein (1204-1503AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transient receptor potential cation channel subfamily M member 2 (Estrogen-responsive element-associated gene 1 protein) (Long transient receptor potential channel 2) (LTrpC-2) (LTrpC2) (Transient receptor potential channel 7) (TrpC7) (Transient receptor potential melastatin 2), TRPM2, EREG1 KNP3 LTRPC2 TRPC7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRPM2

Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA025084ESR1HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.