




CD1A Antibody

Product Code	CSB-PA004889ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06126
Immunogen	Recombinant Human T-cell surface glycoprotein CD1a protein (21-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD1a (T-cell surface antigen T6/Leu-6) (hTa1 thymocyte antigen) (CD antigen CD1a), CD1A
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD1A

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA004889ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.