




# IFNGR1 Antibody

<b>Product Code</b>	CSB-PA011051DSR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P15260
<b>Immunogen</b>	Recombinant Human Interferon gamma receptor 1 protein (270-489AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:200-1:1000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	0.03% Proclin 300,50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Immunology
<b>Target Names</b>	IFNGR1

## Image



### Western blot

All lanes: IFNGR1 antibody at 2µg/ml + Mouse lung tissue  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 55, 22 kDa  
 Observed band size: 55 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.