




# N6AMT1 Antibody

<b>Product Code</b>	CSB-PA010286ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y5N5
<b>Immunogen</b>	Recombinant Human HemK methyltransferase family member 2 protein (1-186AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	HemK methyltransferase family member 2 (EC 2.1.1.-) (M.HsaHemK2P) (N(6)-adenine-specific DNA methyltransferase 1), N6AMT1, C21orf127 HEMK2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	N6AMT1

## Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA010286ESR1HU at dilution of 1:100