





SYTL4 Antibody

Product Code	CSB-PA846595ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C24
Immunogen	Recombinant Human Synaptotagmin-like protein 4 protein (50-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptotagmin-like protein 4 (Exophilin-2) (Granuphilin), SYTL4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SYTL4

Image



Western blot
 All lanes: SYTL4 antibody at 1.36µg/ml + HUVEC whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 77, 40 kDa
 Observed band size: 77 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA846595ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.